

## Product datasheet for **HK207683**

### TUSC3 Human qPCR Template Standard (NM\_006765)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_006765
Synonyms:	D8S1992; M33; MagT2; MRT7; MRT22; N33; OST3A; SLC58A2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209643) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)