

Product datasheet for **HK207707**

TRP1 (TYRP1) Human qPCR Template Standard (NM_000550)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_000550 |
| Synonyms: | b-PROTEIN; CAS2; CATB; GP75; OCA3; TRP; TRP1; TYRP |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP200522) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)