

## Product datasheet for **HK207769**

### UBL5 Human qPCR Template Standard (NM\_024292)

#### Product data:

|                  |  |
|------------------|--|
| Product Type:    | qPCR Template Standards  |
| ACCN:            | NM_024292  |
| Synonyms:        | HUB1   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT   |
| Storage:         | Store at -20°C.  |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.  |



[View online »](#)