

Product datasheet for **HK207784**

PGP9.5 (UCHL1) Human qPCR Template Standard (NM_004181)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004181
Synonyms:	HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207535) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)