

Product datasheet for **HK207794**

UGT (UGT1A9) Human qPCR Template Standard (NM_021027)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021027
Synonyms:	HLUGP4; LUGP4; UDPGT; UDPGT 1-9; UGT-1I; UGT1-09; UGT1-9; UGT1.9; UGT1A9S; UGT1AI; UGT1I
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214005) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)