

Product datasheet for **HK207938**

WDR45 Human qPCR Template Standard (NM_007075)

Product data:

| | |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_007075 |
| Synonyms: | JM5; NBIA4; NBIA5; WDRX1; WIPI-4; WIPI4 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP227408) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)