

Product datasheet for **HK208260**

Cytochrome P450 2D6 (CYP2D6) Human qPCR Template Standard (NM_000106)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000106
Synonyms:	CPD6; CYP2D; CYP2D7AP; CYP2D7BP; CYP2D7P2; CYP2D8P2; CYP2DL1; CYP1D6; P450-DB1; P450C2D; P450DB1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200091) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)