

Product datasheet for **HK208319**

Thyroid Hormone Receptor beta (THRB) Human qPCR Template Standard (NM_000461)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000461
Synonyms:	C-ERBA-2; C-ERBA-BETA; ERBA2; GRTH; NR1A2; PPTH; THR1; THRB1; THRB2; TRbeta
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200441) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)