

Product datasheet for **HK208337**

HNF1 alpha (HNF1A) Human qPCR Template Standard (NM_000545)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000545
Synonyms:	HNF-1A; HNF1; HNF1alpha; HNF4A; IDDM20; LFB1; MODY3; TCF-1; TCF1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200517) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)