

## Product datasheet for **HK208343**

### CD16b (FCGR3B) Human qPCR Template Standard (NM\_000570)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_000570
Synonyms:	CD16; CD16A; CD16b; FCG3; FCGR3; FCGR3A; FCR-10; FCRIII; FCRIIIb
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP200539) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)