

Product datasheet for **HK208486**

NT5C3 (NT5C3A) Human qPCR Template Standard (NM_001002009)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_001002009
Synonyms:	cN-III; hUMP1; NT5C3; P5'N-1; P5N-1; p36; PN-I; POMP; PSN1; UMPH; UMPH1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226396) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)