

## Product datasheet for **HK208775**

### **KREMEN1 Human qPCR Template Standard (NM\_001039571)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_001039571
<b>Synonyms:</b>	FLJ31863; KREMEM1; KREMEN; kringle-coding gene marking the eye and the nose; kringle-containing transmembrane protein 1; kringle containing transmembrane protein 1; KRM1; OTTHUMP00000028977
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP229826) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)