

Product datasheet for **HK208932**

FAAH1 (FAAH) Human qPCR Template Standard (NM_001441)

Product data:

| | |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_001441 |
| Synonyms: | FAAH-1; PSAB |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP205320) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)