

Product datasheet for **HK209223**

SCNN1D Human qPCR Template Standard (NM_002978)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_002978
Synonyms:	dNaCh; ENaCd; ENaCdelta; MGC149710; MGC149711; OTTHUMP0000002084; SCNED; sodium channel, nonvoltage-gated 1, delta; sodium channel, voltage-gated, type I, delta polypeptide
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206584) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)