

Product datasheet for **HK209231**

SFRP4 Human qPCR Template Standard (NM_003014)

Product data:

Product Type: qPCR Template Standards

ACCN: NM_003014

Synonyms: FRP-4; FRPHE; FRZB-2; PYL; sFRP-4

Components: 1 vial of lyophilized SybGREEN qPCR primer (HP206614) mix (1 nmol each primer, sufficient for 200 reactions)
1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA
1 vial containing 1 mL 10x PCR template standard dilution buffer

Quality Control: The kit has been tested to generate satisfactory qPCR data on ABI 7900HT

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)