

Product datasheet for **HK209254**

S6K1 (RPS6KB1) Human qPCR Template Standard (NM_003161)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003161
Synonyms:	p70 S6KA; p70(S6K)-alpha; p70-alpha; p70-S6K; PS6K; S6K; S6K-beta-1; S6K1; STK14A
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP206736) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)