

## Product datasheet for **HK209328**

### H2BC4 Human qPCR Template Standard (NM\_003526)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003526
Synonyms:	dj221C16.3; H2B.1; H2B/l; H2BC6; H2BC7; H2BC8; H2BC10; H2BFL; HIST1H2BC
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207030) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)