

## Product datasheet for **HK209331**

### **H3FH (HIST1H3G) Human qPCR Template Standard (NM\_003534)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | NM_003534   |
| <b>Synonyms:</b>        | H3/h; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C10; H3C11; H3C12; H3FH; HIST1H3G   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (HP207038) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| <b>Storage:</b>         | Store at -20°C.   |
| <b>Stability:</b>       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)