

Product datasheet for **HK209333**

H4C14 Human qPCR Template Standard (NM_003548)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_003548
Synonyms:	FO108; H4; H4-16; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C15; H4F2; H4FN; HIST2H4; HIST2H4A
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207052) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)