

## Product datasheet for **HK209457**

### TRBP (TARBP2) Human qPCR Template Standard (NM\_004178)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004178
Synonyms:	LOQS; TRBP; TRBP1; TRBP2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP229700) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)