

Product datasheet for **HK209491**

RUNX2 Human qPCR Template Standard (NM_004348)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004348
Synonyms:	AML3; CBFA1; CCD; CCD1; OSF-2; OSF2; PEA2aA; PEBP2A1; PEBP2A2; PEBP2aA; PEBP2aA1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225916) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)