

Product datasheet for **HK209564**

eRF1 (ETF1) Human qPCR Template Standard (NM_004730)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004730
Synonyms:	D5S1995; ERF; ERF1; RF1; SUP45L1; TB3-1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208001) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)