

Product datasheet for **HK209573**

Adiponectin (ADIPOQ) Human qPCR Template Standard (NM_004797)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_004797
Synonyms:	ACDC; ACRP30; ADIPQTL1; ADPN; APM-1; APM1; GBP28
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208060) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)