

Product datasheet for **HK209743**

ABI2 Human qPCR Template Standard (NM_005759)

Product data:

| | |
|------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_005759 |
| Synonyms: | ABI-2; ABI2B; AbIBP3; AIP-1; AIP1; argBP1; argBPIA; argBPIB; SSH3BP2 |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP208838) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)