

Product datasheet for **HK210020**

GPR161 Human qPCR Template Standard (NM_007369)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_007369
Synonyms:	FLJ33952; G-protein coupled receptor; G protein-coupled receptor 161; OTTHUMP00000032585; OTTHUMP00000032586; OTTHUMP00000032587; RE2; RE2, FLJ33952
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP227193) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)