

Product datasheet for **HK210026**

Peroxiredoxin 5 (PRDX5) Human qPCR Template Standard (NM_012094)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_012094
Synonyms:	ACR1; AOEB166; B166; HEL-S-55; PLP; PMP20; PRDX6; prx-V; PRXV; SBBI10
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210097) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)