

## Product datasheet for **HK210041**

### **FOXI1 Human qPCR Template Standard (NM\_012188)**

#### **Product data:**

|                  |                                                                                                                                                                                                                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | qPCR Template Standards                                                                                                                                                                                                                                                                       |
| ACCN:            | NM_012188                                                                                                                                                                                                                                                                                     |
| Synonyms:        | FKH10; FKHL10; FREAC-6; FREAC6; HFH-3; HFH3                                                                                                                                                                                                                                                   |
| Components:      | 1 vial of lyophilized SybGREEN qPCR primer (HP210175) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT                                                                                                                                                                                                                      |
| Storage:         | Store at -20°C.                                                                                                                                                                                                                                                                               |
| Stability:       | The primer mix and template standard are stable for one year from date of shipping.                                                                                                                                                                                                           |



[View online »](#)