

Product datasheet for **HK210095**

CD24 Human qPCR Template Standard (NM_013230)

Product data:

| | |
|------------------|--|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_013230 |
| Synonyms: | CD24A |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)