

Product datasheet for **HK210178**

HABP4 Human qPCR Template Standard (NM_014282)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_014282
Synonyms:	IHABP-4; IHABP4; Ki-1/57; SERBP1L
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210761) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)