

Product datasheet for **HK210204**

C7orf54 Human qPCR Template Standard (NM_014411)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_014411
Synonyms:	brain and nasopharyngeal carcinoma susceptibility protein; chromosome 7 open reading frame 54; MGC138346; NAG8; NAG8, NSG-X, MGC138346; nasopharyngeal carcinoma associated gene protein-8; NSG-X
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210866) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)