

Product datasheet for **HK210310**

ICOS Ligand (ICOSLG) Human qPCR Template Standard (NM_015259)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_015259
Synonyms:	B7-H2; B7h; B7H2; B7RP-1; B7RP1; CD275; GL50; ICOS-L; ICOSL; LICOS
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211500) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)