

Product datasheet for **HK210763**

Leukotriene B4 Receptor 2 (LTB4R2) Human qPCR Template Standard (NM_019839)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_019839
Synonyms:	BLT2; BLTR2; JULF2; KPG_004; LTB4-R 2; LTB4-R2; NOP9
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213457) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)