

Product datasheet for **HK210869**

HIST1H2AK (HIST1H2AG) Human qPCR Template Standard (NM_021064)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021064
Synonyms:	H2A.1b; H2A/p; H2AC13; H2AC15; H2AC16; H2AC17; H2AFP; H2AG; HIST1H2AG; pH2A/f
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214024) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)