

Product datasheet for **HK210935**

H4C11 Human qPCR Template Standard (NM_021968)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021968
Synonyms:	dj160A22.2; H4-16; H4/e; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C12; H4C13; H4C14; H4C15; H4F2iv; H4FE; HIST1H4J
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214286) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)