

Product datasheet for **HK210938**

ZNF277 Human qPCR Template Standard (NM_021994)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_021994
Synonyms:	NRIF4; ZNF277P
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214300) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50×10^7 copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)