

Product datasheet for **HK210957**

GI24 (C10orf54) Human qPCR Template Standard (NM_022153)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_022153
Synonyms:	B7-H5; B7H5; C10orf54; DD1alpha; Dies1; GI24; PD-1H; PP2135; SISP1; VISTA
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP214399) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)