

Product datasheet for **HK211203**

KHDC1 Human qPCR Template Standard (NM_030568)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_030568
Synonyms:	bA257K9.4; C6orf147; C6orf148; Em:AC019205.8; KHDC1L; NDG1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP215336) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)