

Product datasheet for **HK211424**

POLR2J2 Human qPCR Template Standard (NM_032958)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_032958
Synonyms:	DNA directed RNA polymerase II polypeptide J-related; HRPB11B; MGC54043; MGC105050; polymerase (RNA) II (DNA directed) polypeptide J2; RPB11, HRPB11B, MGC54043, MGC105050; RPB11b1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP216267) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)