

Product datasheet for **HK211426**

Macrophage inflammatory protein 5 (CCL15) Human qPCR Template Standard (NM_032965)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_032965
Synonyms:	HCC-2; HMRP-2B; LKN-1; LKN1; MIP-1 delta; MIP-1D; MIP-5; MRP-2B; NCC-3; NCC3; SCYA15; SCYL3; SY15
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232796) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)