

## Product datasheet for **HK211512**

### **LY108 (SLAMF6) Human qPCR Template Standard (NM\_052931)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_052931
<b>Synonyms:</b>	CD352; KALI; KAL1b; Ly108; NTB-A; NTBA; SF2000
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP216533) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)