

Product datasheet for **HK211543**

C21orf66 Human qPCR Template Standard (NM_058191)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_058191
Synonyms:	BM-020; chromosome 21 open reading frame 66; FLJ90561; fSAP105; functional spliceosome-associated protein 105; GC-rich sequence DNA-binding factor candidate; GCFC
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP234238) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)