

## Product datasheet for **HK211551**

### **p21 (CDKN1A) Human qPCR Template Standard (NM\_078467)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	NM_078467
<b>Synonyms:</b>	CAP20; CDKN1; CIP1; MDA-6; P21; p21CIP1; SDI1; WAF1
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP200369) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)