

## Product datasheet for **HK211580**

### **TBX1 Human qPCR Template Standard (NM\_080647)**

#### **Product data:**

Product Type:	qPCR Template Standards
ACCN:	NM_080647
Synonyms:	CAFS; CATCH22; CTHM; DGCR; DGS; DORV; TBX1C; TGA; VCF; VCFS
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226193) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)