

Product datasheet for **HK211713**

JNK2 (MAPK9) Human qPCR Template Standard (NM_139069)

Product data:

| | |
|-------------------------|---|
| Product Type: | qPCR Template Standards |
| ACCN: | NM_139069 |
| Synonyms: | JNK-55; JNK2; JNK2A; JNK2ALPHA; JNK2B; JNK2BETA; p54a; p54aSAPK; PRKM9; SAPK; SAPK1a |
| Components: | 1 vial of lyophilized SybGREEN qPCR primer (HP206385) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer |
| Quality Control: | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT |
| Storage: | Store at -20°C. |
| Stability: | The primer mix and template standard are stable for one year from date of shipping. |



[View online »](#)