

## Product datasheet for **HK211762**

### **Kallikrein 11 (KLK11) Human qPCR Template Standard (NM\_144947)**

#### **Product data:**

Product Type:	qPCR Template Standards
ACCN:	NM_144947
Synonyms:	PRSS20; TLSP
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226013) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)