

## Product datasheet for **HK211784**

### MEK7 (MAP2K7) Human qPCR Template Standard (NM\_145185)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_145185
Synonyms:	JNKK2; MAPKK7; MEK; MEK 7; MKK7; PRKMK7; SAPKK-4; SAPKK4
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226278) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)