

Product datasheet for **HK211824**

HMGA1 Human qPCR Template Standard (NM_145904)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_145904
Synonyms:	AL023995; high mobility group AT-hook 1; high mobility group protein I; HMG-I(Y); HMG-R, HMGIY, HMGA1A, MGC4242, MGC4854, MGC12816; Hmga1a; Hmga1b; Hmgi; HMGI(Y); Hmgiy; HMGY; MGC102580
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP232304) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)