

Product datasheet for **HK211859**

SLC39A13 Human qPCR Template Standard (NM_152264)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_152264
Synonyms:	EDSSPD3; LZT-Hs9; SCDEDS; ZIP13
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP217564) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)