

## Product datasheet for **HK212093**

### IL36RN Human qPCR Template Standard (NM\_173170)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_173170
Synonyms:	FIL1; FIL1(DELTA); FIL1D; IL-36Ra; IL1F5; IL1HY1; IL1L1; IL1RP3; IL36RA; PSORP; PSORS14
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP210239) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)