

Product datasheet for **HK212205**

CSPS (SULTIA3) Human qPCR Template Standard (NM_177552)

Product data:

Product Type:	qPCR Template Standards
ACCN:	NM_177552
Synonyms:	HAST; HAST3; M-PST; STIA3; STIA3/STIA4; STIA4; STIA5; STM; TL-PST
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP225807) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.

